



## **Specifications:**

**Stature:** Thigh or Tall Adult Patients

**Range:**

39 to 62 cm

Bariatric (62cm to 124cm) when cascaded

**Bladder Width:** 17cm

**Bladder Volume:**

Restricted =

Non-Restricted=

**Material:**

Medical Grade PVC

Nylon

Latex Free



## **Regulatory:**

AHA, AAMI Sp-10 and IEC-60601 compliant

FDA 510K

## **Coupling:**

Exclusive dual molded tubes form the female connector for use with TEC-Connect system or industry standard Rectus male connectors.

## **Minimum Service Life:**

Inflate 10,000 cycles @ 350mmHg in a restricted mode.

AAMI Sp-10 for reusable = (10,000 @ 330mmHg)

Inflate 100 cycles @ 350mmHg in a non-restrictive mode.

No AAMI specification

Functional hook and loop after 1,000 cycles @ 350mmHg.

AAMI Sp-10 for reusable = (1,000 @ no abuse pressure)

BP cuff transducer is fully functional after 10 washing cycles using automatic machines.